


INSTRUCTIONS FOR USE



Coring
Corneal Ring Segments




Device Description:

CORING - The Corneal Ring Segment Implants are manufactured from Clinical Grade Polymethyl methacrylate (PMMA) and are available from 0.150 mm to 0.400 mm thickness, in 0.050mm increments. The Corneal Ring segments have an arc length from 90 degree to 350 degree. They are lathe cut and tumble polished.

PRODUCT TECHNICAL SPECIFICATION:

Model	Arc Angle	Thickness μm	Width mm	Drawing	Other features
CR-Deg-Th	90°-350°	150-450	0.60		Variable Arc length Thickness
CR6- Deg-Th	90°-350°	150-450	0.60		Variable Arc length Thickness



INDICATIONS:

CORING - Corneal segment Implants are an ophthalmic medical device designed for the reduction or elimination of myopia and astigmatism in patients with keratoconus (KC) so that their functional vision may be restored and the need for a corneal transplant procedure can potentially be deferred.

INTENDED USER:

Ophthalmic surgeons only.

INTENDED PURPOSE:

CORING, Corneal implants are intended for the reduction or elimination of Myopia and Astigmatism in patients with keratoconus, who are no longer able to achieve adequate vision with their contact lenses or spectacle's.

PATIENT TARGET GROUP:

Adult patients of age 18 years and above

INSTRUCTIONS FOR USE

MODE OF ACTION:

When placed in the corneal stroma, outside of the patient's central optical zone, the product reduces the cone by flattening the cornea and for non-central keratoconus, repositions the cone centrally. Corneal segments are designed to be placed in the periphery of the cornea, at approximately two-thirds depth, and are surgically inserted through a small radial incision in the corneal stroma. The placement of the incision will be typically temporal at the axis of positive cylinder; however, may vary depending on the location of the cone and the amount of keratoconus present in the specific eye to be treated. The Corneal segments are to be placed equidistant on each side of the incision.

To reduce the myopia and the irregular astigmatism induced by keratoconus, two Corneal segments may be implanted depending on the patient's preoperative manifest refraction spherical equivalent (MRSE), the orientation of the cone and the degree of asymmetric astigmatism. The product is designed with a fixed outer diameter and width. The Corneal segments have two positioning holes, located at each end of the segment, to aid in surgical manipulation. CORING segments are available for Optic zone diameter 5 mm and 6 mm.

II. Typical Treatment Nomogram (see diagram below).

The Corneal segments treatment nomogram for keratoconus is based on the clinical results from implanting the 0.350 mm, 0.400 mm and 0.450 mm thicknesses of segments in keratoconus patients. The Lower thickness is proposed to be used to treat those keratoconus patients who require only a minimal amount of corneal flattening and/or corneal stabilization, who are contact lens intolerant and who can no longer be effectively corrected with spectacles to restore their functional vision. The specific determination of which thicknesses of the corneal segments to implant is dependent upon the nomogram and several variables: the most significant being the patient's preoperative manifest refraction spherical equivalent, the location of the cone and the degree of asymmetric astigmatism present.

- The surgical technique for keratoconus is similar to the standard corneal segment surgical technique used for low myopia, except that the location of the incision is often placed temporally.
- Pachymetry is to be measured during surgery at the peripheral location of the entry incision. All patients shall have a corneal thickness of 450 microns or greater at the proposed incision site. The incision depth should be at 68% of the corneal thickness measured at the peripheral location of the entry incision.
- A temporal incision is typically used depending on the location of the astigmatism. Some corneal surgeons place the entry incision in the same meridian as the axis of positive cylinder.
- Note: Patients may be allowed to manually adjust the axis of cylinder at the Phoropter to achieve the clearest subjective image. This technique may work better than the Jackson cross-cylinder technique for patients with irregular astigmatism.
- An asymmetrical cone is treated with a thinner Corneal segment placed superiorly and a thicker Corneal segment placed inferiorly.
- Global and central cones are typically treated using two segments of the same thickness.
- The thicker Corneal segment is typically placed to correspond to the keratoconus cone (inferiorly) to lift the cone and produce the maximum flattening effect and the thinner Corneal segment is placed in the opposite corneal half (superiorly) to counterbalance the thicker segment and to flatten the rest of the corneal surface across the visual axis.

There are two primary criteria used in determining the surgical nomogram related to the use of Corneal segments Implants for keratoconus. The first criterion is whether the cone is centered or decentered (asymmetric cone). Keratoconus, which is centrally present as determined by a topographical map, will require two corneal segments of the same thickness.

The thickness of the Corneal segments to be used is determined based on the preoperative spherical equivalent of less than or equal to -3.00 D or greater than -3.00 D. This would also apply to global keratoconus in which the keratoconus is central, but its circumference extends beyond 5.0 mm from the center.

INSTRUCTIONS FOR USE

As for asymmetric cones, it is necessary to evaluate the degree to which the cone is decentered. This is done by reviewing a topographical map of the cornea. Moderate asymmetry exists when the cone is off-center in placement at the 3.0 mm ring on the topographical map. High asymmetric cones are typically 5.0 mm or more off-center, as exhibited on a topographical map. In each of these cases, two different thicknesses used will depend on whether the preoperative spherical equivalent is less than or equal to -3.00 D or greater than -3.00 D. The thicker Corneal segment is recommended to be placed inferiorly and the thinner Corneal segment is to be recommended be placed superiorly.




The recommended corneal segment placement and thickness nomogram to be used for keratoconus is presented below:

Recommended Corneal segment Placement Nomogram for Keratoconus	
Type of Keratoconus	Product Configuration
Asymmetrical Cone	Two Corneal Segments of different thicknesses implanted (One thinner Corneal segment placed superiorly – One thicker Corneal segment placed inferiorly)
Global Cone	Two corneal segments of the same thickness implanted (One Corneal segment Placed superiorly – One Corneal segment placed inferiorly)
Central Cone	Two corneal segments of the same thickness implanted (One Corneal segment Placed superiorly – One Corneal segment placed inferiorly)

Recommended corneal segment Thickness Nomogram for keratoconus		
Type of cone	Preop MRSE < 3.00 D	Preop MRSE > 3.00 D
Asymmetrical Cone: Moderate asymmetry High asymmetry	0.300 mm 0.350 mm/0.400 mm 0.400 mm	0.350 mm 0.400 mm/0.450 mm 0.450 mm
Global cone	0.400 mm / 0.400 mm	0.450 mm/0.450 mm
Central cone	0.400 mm / 0.400 mm	0.450 mm/0.450 mm

Some of the different product configurations recommended for the treatment of keratoconus are illustrated below in Diagrams 1- 3:

INSTRUCTIONS FOR USE

<p>Diagram 1 Asymmetrical Cone-</p> 	<p>For patients with a preoperative MRSE \leq 3.00 D, one thinner Corneal segment (e.g., 0.300 mm) is placed superiorly and one thicker Corneal segment is placed inferiorly (e.g., 0.350mm), temporal incision placement is used (9:00 for right eye and 3:00 for left eye).</p>
<p>Diagram 2 Asymmetrical Cone</p> 	<p>For patients with a preoperative MRSE $>$ 3.00 D, one Corneal segment (e.g., 0.300 mm) is placed superiorly and one thicker corneal segment is placed inferiorly. (e.g., temporal incision used (9:00 for right eye and 3:00 for left eye)</p>
<p>Diagram 3 Global Cone/Central Cone</p> 	<p>For patients with a preoperative MRSE $>$ 3.00D, two segments of the same thickness are placed – one segment is placed superiorly and one segment is placed inferiorly, temporal incision placement is used (9:00 for right eye and 3:00 for left eye).</p>

Indications:

Corneal segments (Implants) are intended for the reduction or elimination of myopia and Astigmatism in patients with keratoconus, who are no longer able to achieve adequate vision with their contact lenses or spectacles, so that their functional vision may be restored and the need for a corneal transplant procedure may potentially be deferred.

The specific subset of keratoconus patients proposed to be treated with Corneal segment Implants are those patients:

- Who have experienced a progressive deterioration in their vision, such that they can no longer achieve adequate functional vision on a daily basis with their contact lenses or spectacles.
- Who are 21 years of age or older.
- Who have clear central corneas.
- who have a corneal thickness of 450 microns or greater at the proposed incision site; and
- who have corneal transplantation as the only remaining option to improve their functional vision.

CLINICAL BENEFITS:

- The clinical benefit of the implantation of Corneal Ring Segment is to avoid or postpone Corneal Transplant.
- CORING Segments improves patient's visual acuity.
- Avoid spectacle dependence and use of uncomfortable hard contact lens.

SIDE EFFECTS:

- Foreign Body Sensation
- Corneal Thinning

RESIDUAL RISK:

Some patients may experience loss of contrast sensitivity post operatively.

Contraindications:

Corneal segment (implants) for keratoconus is Contraindicated:

INSTRUCTIONS FOR USE

- In patients who have abnormally thin corneas or who have a corneal thickness of 499 microns or less at the proposed incision site.
- In patients with collagen vascular, autoimmune or immunodeficiency diseases.
- In pregnant or nursing women.
- In the presence of ocular conditions, such as recurrent corneal erosion syndrome or corneal dystrophy, that may predispose the patient to future complications; or
- In patients who are taking one or more of the following medications: Isotretinoin (Accutane¹); amiodarone hydrochloride (cordarone²).
- In patients who have Fuchs' dystrophy or have suspect endothelial cell deficiency are not candidates for corneal segment Implants.

Warnings:

- Some patients with large, dilated pupil diameters (≥ 7.0 mm) are predisposed to low light visual symptoms postoperatively and should be appropriately advised.
- The long-term effect of corneal segment Implants on endothelial cell density has not been established. Central endothelial cell density loss for myopic eyes implanted with 0.350 mm corneal segment implants was $1.4\% \pm 3.9\%$.
- Patients who have the 0.400 mm or 0.450 mm corneal segment Implants should be monitored on a more frequent basis in order to detect peripheral endothelial cell loss and/or related corneal edema. These patients should be advised about the potential risk of developing peripheral endothelial cell loss and/or corneal edema, possibly requiring corneal transplantation, with potential unforeseen impact of peripheral endothelial cell loss on the success of future PKP.

Annualizes Percent Change in Endothelial Cell density, 0.400 mm, and 0.450 mm Corneal segment Implants

Region	Central	6.00 Peripheral	10.00 Peripheral
N	17	12	15
Mean \pm SD (%/year)	-1.63% \pm .79%	-2.22% \pm 0.86%	-1.85% \pm 0.56%
95% Confidence Interval (%/year)	-1.99%-1.26%	-2.71%: -1.73%	-2.13%: -1.56%

- Under mesopic conditions, patients may experience some loss in contrast sensitivity at low spatial frequencies (1.5 cycle per degree).

Precautions

Use of the vacuum Centring Guide subjects the eye to increased intraocular pressure. Continuous application of vacuum should be limited to 3 minutes or less and to no more than 750 m Bar. If it is necessary to reapply the Vacuum Centring Guide, wait 5 minutes to allow normal vascular perfusion of the eye to occur before re-establishing suction.

- Corneal segment Implants are not recommended in patients with systemic diseases likely to affect wound healing, such as insulin-dependent diabetes or severe atopic disease.
- It is critical that the corneal segment Implants are properly oriented during the surgical procedure so that the corneal segment are not left beneath the incision area in order to avoid potential wound healing issues. The alignment of the implanted corneal segments is to be verified against the placement marks made on the cornea at the beginning of the corneal segment implantation procedure. If the corneal segments are not properly oriented operatively, the segments should be immediately re-positioned.
- Some patients (7/334 or 2.1%) who received the 0.400 mm and 0.450 mm Corneal segment Implants for the treatment of myopia experienced a temporary decrease of two or more lines of BSCVA during

INSTRUCTIONS FOR USE

the U.S. myopia trial. The BSCVA was 20/20 or better upon resolution for all patients and the BSVCA for all patients returned to within 2-8 letters of their preoperative BSCVA

- It is recommended that the corneal thickness at the 6 mm to 9 mm optical zone be at least 350 microns to allow for an adequate amount of corneal tissue above the Corneal segment Implants.
- Corneal segment Implants are not recommended in patients with a history of ophthalmic herpes simplex or Herpes zoster.
- Corneal segment Implants are not recommended in patients who are taking sumatriptan (Imitrex³) for migraine headaches.
- A temporary decrease in central corneal sensation has been noted in some patients.
- The safety and effectiveness of alternative refractive procedures or a corneal transplant procedure following the removal of Corneal segment Implants have not been established.
- Corneal segment Implants are intended for single use only; do not reuse or resterilize. In the event that different thicknesses of corneal segment Implants are used during a procedure, please discard the unused segments.
- The safety of corneal segment Implants for keratoconus have NOT been established:
 - In patients with progressive myopia or astigmatism, nuclear sclerosis or other crystalline lens opacity, corneal abnormality, or previous corneal surgery or trauma.
 - For patients under 21 years of age.
 - For corneas with a central thickness less than 480 microns, or peripheral thickness less than 570 microns; or
 - In long-term use.

Patient Instructions and Identification Card patient Instructions

- If patients wear contact lenses, they should be instructed to stop wearing them 2-3 weeks before their preoperative examination in order to obtain an accurate refraction.
- Patients should be instructed on the importance of using all medications as directed.
- Patients should be instructed to use the night time eye shield as directed to avoid injuring their surgery eye during sleep.
- Patients should be instructed not to rub their surgery eye for the first six months after the procedure. This is important to promote proper healing of the incision.
- Patients should be instructed to avoid getting tap water in their surgery eye for the first few weeks following the procedure.
- Patients should be advised to avoid swimming in pools for the first week after the procedure. Lake and ocean swimming should be avoided for the first month. Protective goggles should be worn at all times when swimming
- Patients should be instructed to contact you immediately if they experience any pain, discomfort, have a sensation that something is in their eye or experience a change in their vision after the initial postoperative recovery period (typically 7 days). Patients should be instructed to report any unusual symptoms that could be associated with prolonged topical steroid use, if applicable.

Identification Card

A Patient Identification Card is enclosed in the CORING Implants product package.

The implant card is to be completed by the healthcare provider. A product traceability label also supplied with this device must be affixed to the implant card after completion of the surgery. Additional labels can be used for the patient file or clinical follow-up. Completed implant cards must be provided to the patient post-procedure. Patients should be instructed to keep the card as a permanent record of their implant. The patient should also be instructed.

to show the card to any eye practitioner he or she may see in the future.

Conformance to Standards

CORING - Corneal segment Implants have been designed, manufactured, and distributed in conformance with requirements of the Quality System management ISO 13485:2016 and other relevant ISO standards, and the Medical Device Regulation (MDR) 2017/745.

Medical Device Reporting

Doc. No : EXL/IFU/06
Rev.No MDR : 00
Rev.Date : 01/2024

INSTRUCTIONS FOR USE

Users should report the serious incident with medical device information to the manufacturer and/or to the national competent authority depending on the national practice. Once corrective (or other) action is identified from the Manufacturer, hospital administrators, medical practitioners, and other health-care professionals, and USER representatives responsible for the maintenance and the safety of MEDICAL DEVICES, can take the necessary steps. Such steps should, where practicable, be taken in co-operation with the MANUFACTURER. For the purposes of Medical Devices Vigilance System in member states are represented by appointed National Competent Authorities, their vigilance contacts points being listed on the European Commission web site: http://ec.europa.eu/growth/sectors/medicaldevices/contacts/index_en.html

How Supplied

The CORING, Corneal segments Implant are supplied in a plastic case contained within a heat-sealed Tyvek peel pouch and terminally sterilized using Ethylene Oxide (EO). The contents of the pouch are sterile unless the package is opened or damaged.

Return Goods Policy

For information on returning any damaged product, contact your local representative or call Excel optics Pvt Ltd for return authorization and full policy information. All products returned to Excel optics Pvt Ltd must be accompanied by a Return Goods Authorization Number.

LIMITATION OF WARRANTY AND LIABILITY

EXCEL OPTICS (P) LIMITED accepts no liability for any injury suffered to patients as a result of: Any implantation method or technique used by a physician to implant the CORING Corneal Ring Implant.

Any prescription, and use of the CORING for any Individual patient or patient's conditions.

Excel Optics (P) Limited, makes no expressed or implied.

Warranties in connection with the sale of the CORING – Corneal Ring Segments.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE:

Summary of Safety and Clinical Performance (SSCP) of this device will be available on the EUDAMED website <https://ec.europa.eu/tools/eudamed> after notified body acceptance of SSCP.

Note: Because the Corneal Ring and the packaging materials are plastic, the Corneal Ring may pick up an electrostatic charge upon opening the package. The Corneal Ring should be carefully examined to ensure that particles have not been attracted to it.

Implantation of Corneal Ring should not be performed in patients under the age of 18.

CAUTION:

GOVERNMENT LAWS RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

DISPOSAL

Discarded Corneal Ring (used or unused) are classified as medical (clinical) waste that harbours a potential infection or microbial hazard and must be disposed of accordingly.

Do not dispose of damaged or explanted device or its packing with household trash. Disposal of devices and their packaging is considered a biohazard.

Follow local regulatory guidelines for disposing of medical devices and their packaging safely.

STORAGE CONDITIONS:

Store between 5°C to 40°C temperature.





















EXPIRATION DATE:

The expiration date on the Corneal Ring package is the sterility expiration date.

Do not use the Corneal Ring after the expiration date.

Symbols and their Explanations

INSTRUCTIONS FOR USE

	Authorized Representative in the European community		Do Not re-use
	Sterilized using Ethylene Oxide		Do not re-sterilize
	Storage Condition between 5°C to 40°C		Keep away from sunlight
	Manufacturer		Consult Instructions for Use
	Sterile Batch No		Keep Dry
	Unique device identifier		Medical Device
	Serial Number		Single sterile barrier system with protective packaging inside
	Model Number		Expiry Date
	Do not Use If package is Damaged		Date of Manufacturer
	Country of Manufacture		
	Ultraviolet		

	<p>MED DEVICES LIFESCIENCES B.V. Keizersgracht 482, 1017EG Amsterdam, The Netherlands Phone: +31 – 202254558, Email: info@meddevices.net</p>
---	--

<p> EXCEL OPTICS (P) LIMITED New No: 31A & 31B, Old. No: 15A & 15 B, 2nd Street, Logan than Nagar, Choolaimedu, Chennai – 600 094.INDIA E-mail: excel_factory@hotmail.com www.exceloptics.net</p>	<p>CE 1434 ISO 13485:2016 Certified Company</p>
--	--